

## Turnkey cleanroom with isolators

### Blackpool Victoria Hospital

Pharmacy aseptic unit for cytotoxics and CIVAS



#### Principle features of project:

- ❑ Turnkey package with room build by Bassaire Projects and equipment by Envair
- ❑ Total project value: over £400k
- ❑ Close co-operation with client on validation and documentation
- ❑ Full IQ OQ documentation including Factory Acceptance Tests (FAT) and Site Acceptance Tests (SAT)

#### Principle features of facility:

- ❑ Two Envair Pharm-Assist positive pressure isolators for CIVAS
- ❑ Two Envair CDC F negative pressure isolators for cytotoxics
- ❑ Full fresh air / full dump for cytotoxic room to eliminate the risk of cross contamination
- ❑ Controls interface with site BMS system, with dedicated alarming to email system
- ❑ HTM specifications for HVAC, benching, furniture, finishes, etc
- ❑ Future proofing for potential changes to EU GMP grades